

Interton - Essential Hearing

Founded in Germany in 1962, Interton today is part of one of the world's largest hearing instrument organizations, the GN Group. We use their proven core technology in a way that falls in line with our own way of doing things. We believe that most people need a straight-forward, functional and easy-to-use hearing solution at a fair price. So that's what we give them.

Interton is a trademark of GN Hearing A/S. © 2016 GN Hearing A/S. All rights reserved. Apple, the Apple logo, iPhone, iPad, iPod touch and FaceTime are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android and Google Play are trademarks of Google Inc. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories. Bluetooth is a trademark of Bluetooth SIG, Inc. The Interton EasyHearing App is compatible with iPhone, iPad and iPod touch using iOS 7.0 or later. The Android version of Interton EasyHearing App is available on selected Android devices running version 2.3 or later.

ENHANCE YOUR HEARING EXPERIENCE WITH WIRELESS ACCESSORIES

M102567-GB-16.09-Rev.A



Make every hearing experience count

Even the best hearing aids can sometimes use a boost. Restaurants, noisy parties or busy events where you are listening to someone from a distance can be especially challenging. Or when you are driving a car, trying to tune in to conversations from the back seat.

Interton's second generation of industry-leading 2.4 GHz wireless accessories make all of life's listening situations clearer and more pleasurable. They help you pick up on quiet details and subtle nuances so, you can have confidence that you won't miss a thing.

Stylish, easy to use and simple to set up, our wireless accessories connect to your phone, TV, computer or stereo and receive speech and other essential sounds directly in your hearing aids.

NEW Mini Microphones

Transmit speech or music directly to your wireless hearing aids so you can hear every word - even at a distance.

TV Streamer 2

Stream sounds from your TV directly to your hearing aids, turning them into wireless stereo headphones.

Remote Control 2

Adjust the volume in either one or both hearing aids and change listening programs effortlessly.

Phone Clip 2

Clips discreetly to your clothes and connects your hearing aids directly to a Bluetooth® enabled mobile phone.

Small, portable and extremely easy to use

The latest additions to our range of wireless accessories, Interton's Micro Mic and Multi Mic let you hear more of everything. Whether enjoying a family dinner, being in a meeting or listening from a distance, simply ask the speaker to clip on one of our Mini Mics. The Mini Mics stream the speaker's words directly to your wireless hearing aids, so you can hear everything clearly, even when a lot of background noise is present. You can also plug them into your music player to enjoy music anywhere.

Interton Micro Mic

Ideal for non face-to-face conversations, the Micro Mic pairs to your Interton hearing aids in seconds. Clip it on to the speaker's clothing and start listening - even if you are up to 25 metres away in clear line of sight.

Interton Multi Mic

Works like the Micro Mic with the same range of up to 25 metres, but also doubles as a table microphone, so you can hear everyone sitting around it. Connect with loop and FM systems, and use the mini-jack input for streaming audio to your hearing aids.



TV is just better together

Whether you prefer drama, thrillers or comedy, TV is almost always better when you can enjoy it with the ones you love - at a volume level everyone can agree upon. Interton TV Streamer 2 lets you stream sound directly to your wireless hearing aids so everyone can watch TV at a comfortable sound level and talk together at the same time.

The TV Streamer 2 streams superb stereo sound up to seven metres and reconnects automatically if you leave the room and come back within five minutes. It works with any digital or analogue TV or audio source.

Easy overview. Full control

The easiest way to control your hearing aids, TV Streamer 2 or any other wireless accessories is with the Interton Remote Control 2. You can adjust the volume in one or both hearing aids, change listening programmes or easily switch between streaming sources.

With its large LCD display, the Remote Control 2 shows your listening programmes, volume and its battery level at a glance. It has large, intuitive buttons that ensure easy operation. You can balance the hearing aid volume with the volume of your streaming sources and block out sounds with the mute function.



Remote Control 2



TV Streamer 2



Smooth and focused phone calls

You no longer need to struggle and strain when you're on your mobile phone in noisy places. The Interton Phone Clip 2 lets you stream calls directly to your wireless hearing aids so you can enjoy crystal-clear phone conversations - no matter where you are.

With the Phone Clip 2 you can also use the EasyHearing app to change programmes and adjust volume on your hearing aids discreetly and hassle free. And with a pair of Centro 6 hearing aids you can access everything from music to FaceTime calls or GPS instructions and control your hearing aids directly from your iPhone.

Ask your local hearing care professional for more information about Interton's many accessory options.



EasyHearing app



Phone Clip 2

